



Recombinant Rhizobium radiobacter Protein virD3 (virD3)

Product Code	CSB-MP357778RKV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09816
Product Type	Recombinant Protein
Immunogen Species	Rhizobium radiobacter (Agrobacterium tumefaciens) (Agrobacterium radiobacter)
Purity	>85% (SDS-PAGE)
Sequence	MANGQFTIRS AGPASVGLTG ERRGAASASS SALSINVQGDV RDRPIPTSSS RLPNAAILRD SSEIASTGLR YVAATLHWSE IAPLSLINKN DLAPAAIDFE TRNDAGNVTA KVGRAVPVPK QGRLGKMLAS VPLSTRISRV NSDRRLPTDA ENHPETRDPR KGRGSHGVTP ALHEKIGTAF VRRLRKIRTV LFAVAARLGA R
Source	Mammalian cell
Target Names	virD3
Protein Names	Recommended name: Protein virD3
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.