



Recombinant Salmonella typhimurium Cysteine synthase A (cysK)

Product Code	CSB-EP358053SXB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A1E3
Product Type	Recombinant Protein
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Purity	>85% (SDS-PAGE)
Sequence	SKIYEDNSL TIGHTPLVRL NRIGNGRILA KVESRNPSFS VKCRIGANMI WDAEKRGVVK PGVELVEPTS GNTGIALAYV AAARGYKLTL TMPETMSIER RKLLKALGAN LVLTEGAKGM KGAIQKAEI VASDPQKYLL LQQFSNPANP EIHEKTTGPE IWEDTDGQVD VFISGVGTGG TLTGVTYRIK GTKGKTDLIT VAVEPTDSPV IAQALAGEEI KPGPHKIQQI GAGFIPGNLD LKLIDKVVGI TNEEAISTAR RLMEEEGILA GISSGAAVAA ALKLQEDEF TNKNIVVILP SSGERYLSTA LFADLFTEKE LQQ
Source	E.coli
Target Names	cysK
Protein Names	Recommended name: Cysteine synthase A Short name= CSase A EC= 2.5.1.47 Alternative name(s): O-acetylserine (thiol)-lyase A Short name= OAS- TL A O-acetylserine sulfhydrylase A
Expression Region	2-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.