



Recombinant Salmonella typhimurium Cysteine synthase A (cysK)

| | |
|--------------------------|---|
| Product Code | CSB-YP358053SXB |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P0A1E3 |
| Product Type | Recombinant Protein |
| Immunogen Species | Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720) |
| Purity | >85% (SDS-PAGE) |
| Sequence | SKIYEDNSL TIGHTPLVRL NRIGNGRILA KVESRNPSFS VKCRIGANMI WDAEKRGVVK PGVELVEPTS GNTGIALAYV AAARGYKLTL TMPETMSIER RKLLKALGAN LVLTEGAKGM KGAIQKAEI VASDPQKYLL LQQFSNPANP EIHEKTTGPE IWEDTDGQVD VFISGVGTGG TLTGVTYRIK GTKGKTDLIT VAVEPTDSPV IAQALAGEEI KPGPHKIQGI GAGFIPGNLD LKLIDKVVGI TNEEAISTAR RLMEEEGILA GISSGA AVAA ALKLQEDES F TNKNIVVILP SSGERYLSTA LFADLFTEKE LQQ |
| Source | Yeast |
| Target Names | cysK |
| Protein Names | Recommended name: Cysteine synthase A Short name= CSase A EC= 2.5.1.47 Alternative name(s): O-acetylserine (thiol)-lyase A Short name= OAS- TL A O-acetylserine sulfhydrylase A |
| Expression Region | 2-323 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |