



Recombinant Human Interleukin-1 receptor-associated kinase-like 2 (IRAK2)

Product Code	CSB-MP011810HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43187
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MACYIYQLPS WVLDDLCRNM DALSEWDWME FASYVITDLT QLRKIKSMER VQGV SITREL LWWWGMRQAT VQQLVDLLCR LELYRAAQII LNWKPAPPEIR CIPAF PDSV KPEKPLAASV RKAEDEQE EG QPVRMATFPG PGSSPARAHQ PAFLQPPEED APHSLRSDLP TSSDSKDFST SIPKQEKLLS LAGDSLFWSE ADVVQATDDF NQNRKISQGT FADVVRGHRH GKPFVFKKLR ETACSSPGSI ERFFQ AELQI CLRCCHPNVL PVLGFCAARQ FHSFIYPYMA NGSLQDRLQG QGGSDPLPWP QRVSICSGLL CAVEY LHGLE IHSNVKSSN VLLDQNLTPK LAHPMAHLCP VNKRSKYTMM KTHLLRTSAA YLPEDFIRVG QLTRKVDIFS CGIVLAEVLT GIPAMDNNRS PVY LKDLLS DIPSSSTASLC SRKTGVENVM AKEICQKYLE KGAGRLPEDC AEALATAACL CLRRRNTSLQ EVCGSVAAVE ERLRGRETL L PWSGLSEG TG SSSNTPEETD DVDNSSLDAS SSMSVAPWAG AATPLLPTEN GEGRLRVIVG READSSSEAC VGLEPPQDVT ETSWQIEINE AKRKL MENIL LYKEEKVDSI ELFGP
Source	Mammalian cell
Target Names	IRAK2
Protein Names	Recommended name: Interleukin-1 receptor-associated kinase-like 2 Short name= IRAK-2
Expression Region	1-625
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	IRAK2 encodes the interleukin-1 receptor-associated kinase 2, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. IRAK2 is reported to participate in the IL1-induced upregulation of NF-kappaB.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.