



Recombinant Virulence regulon transcriptional activator virF (virF)

Product Code	CSB-MP358278SZB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A2T1
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MMDMGHKNKI DIKVRLHNYI ILYAKRCSMT VSSGNETLTI DEGQIAFIER NIQINVSIIK SDSINPFEII SLDRNLLLSI IRIMEPIYSF QHSYSEEKRG LNKKIFLLSE EEVSIDLFKS IKEMPFGKRK IYSLACLLSA VSDEEALYTS ISIASSLSFS DQIRKIVEKN IEKRWRLSDI SNNLNLSEIA VRKRLESEKL TFQQILLDIR MHHAAKLLL N SQSYINDVSR LIGISSPSYF IRKFNEYGI TPKKFYLYHK KF
Source	Mammalian cell
Target Names	virF
Protein Names	Recommended name: Virulence regulon transcriptional activator virF
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.