



# Recombinant Human Zinc finger and BTB domain-containing protein 43 (ZBTB43)

<b>Product Code</b>	CSB-EP026335HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O43298
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEPGTNSFRV EFPDFSSTIL QKLNQQRQQG QLCDVSIVVQ GHIFRAHKAV LAASSPYFCD QVLLKNSRRI VLPDVMNPRV FENILLSSYT GRLVMPAPEI VSYLTAASFL QMWHVVDKCT EVLEGNPTVL CQKLNHGSDH QSPSSSSYNG LVESFELGSG GHTDFPKAQE LRDGENEEES TKDELSSQLT EHEYLPSSNS TEHDLRSTEM ASQDGEEGAS DSAEFHYTRP MYSKPSIMAH KRWIHVKPER LEQACEGMDV HATYDEHQVT ESINTVQTEH TVQPSGVEED FHIGEKKVEA EFDEQADESN YDEQVDFYGS SMEEFSGERS DGNLIGHRQE AALAAGYSEN IEMVTGIKKEE ASHLGFSATD KLYPCQCGKS FTHKSQRDRH MSMHLGLRPY GCGVCGKKFK MKHHLVGHMK IHTGIKPYEC NICAKRFMWR DSFHRHVTSC TKSYEAAKAE QNTTEAN
<b>Source</b>	E.coli
<b>Target Names</b>	ZBTB43
<b>Protein Names</b>	Recommended name: Zinc finger and BTB domain-containing protein 43 Alternative name(s): Zinc finger and BTB domain-containing protein 22B Zinc finger protein 297B ZnF-x
<b>Expression Region</b>	1-467
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.