



Recombinant Human CBP80/20-dependent translation initiation factor (CTIF)

Product Code	CSB-BP012191HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43310
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MENSSAASAS SEAGSSRSQE IEELERFIDS YVLEYQVQGL LADKTEGDGE SERTQSHISQ WTADCSEPLD SSCSFSRGRA PPQQNGSKDN SLDMLGTDIW AANTFDSFSG ATWDLQPEKL DFTQFHRKVR HTPKQPLPHI DREGCGKGKL EDGDGINLND IEKVLPWQGG YHPMPHEVEI AHTKKLFRRR RNDRRRQQR PGGNKPQQHG DHQPGSAKHND RDHQKSYQGG SAPHPSGRPT HHGYSQNRRW HHGNMKHPPG DKGEAGAHNRN AKETMTIENP KLEDTAGDTG HSSLEAPRSP DTLAPVASER LPPQQSGGPE VETKRKDSIL PERIGERPKEI TLLQSSKDRL RRRLKEKDEV AVETTTTPQQN KMDKLEILN SMRNNSSDVD TKLTTFMEEA QNSTNSEEML GEIVRTIYQK AVSDRSFAFT AAKLCDKMAL FMVEGTFKFRS LLLNMLQKDF TVREELQQQD VERWLGFITF LCEVFGTMRS STGEPFRVLV CPIYTCLREL LQSQDVKEDA VLCCSMELQS TGRLLLEEQLP EMMTELLASA RDKMLCPSES MLTRSLLELV IELHANSWNP LTPPITQYYN RTIQKLTA
Source	Baculovirus
Target Names	CTIF
Protein Names	Recommended name: CBP80/20-dependent translation initiation factor
Expression Region	1-598
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.