




# Recombinant Bordetella bronchiseptica Azurin (BB3856)

<b>Product Code</b>	CSB-YP358313BTY
<b>Abbreviation</b>	Recombinant Bordetella bronchiseptica Azurin protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0A321
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bordetella bronchiseptica (strain ATCC BAA-588 / NCTC 13252 / RB50) (Alcaligenes bronchisepticus)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	AECVSDIAGTDQMQFDKKAIEVSKSCKQFTVNLKHTGKLPNVMGHNWVLT KTADMQAVEKDGIAAGLDNQYLKAGDTRVLAHTKVLGGGESDSVTFDVAKLAA GDDYTFFCSPFGHGALMKGTLKLVD
<b>Research Area</b>	Others
<b>Source</b>	Yeast
<b>Target Names</b>	BB3856
<b>Protein Names</b>	Recommended name: Azurin
<b>Expression Region</b>	22-150aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	15.8 kDa
<b>Protein Length</b>	Full Length of Mature Protein

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

