



Recombinant Human Eukaryotic translation elongation factor 1 epsilon-1 (EEF1E1)

Product Code	CSB-MP007432HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43324
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AAAAELSLLEKSLGLSKGN KYSAQGERQI PVLQTNNGPS LTGLTTIAAH LVKQANKEYL LGSTAEKAI VQQWLEYRVT QVDGHSSKND IHTLLKDLNS YLEDKVYLTG YNFTLADILL YYGLHRFIVD LTVQEKEKYL NVSRWFCHIQ HYPGIRQHLS SVVFIKNRLY TNSH
Source	Mammalian cell
Target Names	EEF1E1
Protein Names	Recommended name: Eukaryotic translation elongation factor 1 epsilon-1 Alternative name(s): Aminoacyl tRNA synthetase complex-interacting multifunctional protein 3 Elongation factor p18 Multisynthase complex auxiliary component p18
Expression Region	2-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a multifunctional protein that localizes both in the cytoplasm and in the nucleus. In the cytoplasm, the encoded protein is an auxiliary component of the macromolecular aminoacyl-tRNA synthase complex. However, its mouse homolog has been shown to translocate to the nucleus in response to DNA damage and plays a positive role in ATM/ATR-mediated p53 activation. Alternative splicing results in multiple transcript variants.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.