



Recombinant Thiol:disulfide interchange protein DsbA-like (c3786)

Product Code	CSB-YP358545EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A4L7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	FTE GTDYMVLEKP IPNADKTLIK VFSYACPCFCY KYDKAVTGPV SEKVKDIVAF TPFHLETKGE YGKQASEVFVA VLINKDKAAG ISLFDANSQF KKAKFAYYAA YHDKKERWSD GKDPAAFIKT GLDAAAGMSQA DFEAALKEPA VQETLEKWKA SYDVAKIQGV PAYVVNGKYL IYTKSIKSID AMADLIRELA SK
Source	Yeast
Target Names	dsbL
Protein Names	Recommended name: Thiol:disulfide interchange protein DsbA-like
Expression Region	28-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.