



Recombinant Rabbit Hepatic triacylglycerol lipase (LIPC)

Product Code	CSB-EP012974RB-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O46559
Product Type	Recombinant Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Purity	>85% (SDS-PAGE)
Sequence	GQSLRPEP FGRRARVTAT KKTLLLETETR FLLFKDKANK GCQIRLHHAD TLQECGFNSS LPLVMIVHGW SVDGLLESWI WQMVAALKSQ PARPVNVGLV DWISLAHSHY AVAVRNARLV GQEVALLQW LEESAPFSRS NVHLIGYSLG AHVAGFAGSY ISGKHKIGRI TGLDAAGPLF EGTSASDRLS PDDATFVDI HTFTREHMGL SVGIKQPVGH YDFYPNGGSF QPGCHFLELY KHIAQHGLNA LSQTIKCAHE RSVHLFIDSL LHPSMQSTAY QCSDMDSFSQ GLCLGCTKGR CNTLGYHIRQ EPLSKGKRLF LVTQAQSPFR VYHYQFKIQF INQIEKPLEP TFTMSLLGTK EEMQKIPITL GEGITSNKTY SFLITLNLDI GELMVIKFKW ENSAVWANVW NTVQTIIPWG IKPRNSGLIL KTIRVKAGET QQRMTFCSEN MDDLQLHPTQ EKNFVRCEVN PKKLKLIK
Source	E.coli
Target Names	LIPC
Protein Names	Recommended name: Hepatic triacylglycerol lipase Short name= HL Short name= Hepatic lipase EC= 3.1.1.3 Alternative name(s): Lipase member C
Expression Region	23-499
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	LIPC encodes hepatic triglyceride lipase, which is expressed in liver. LIPC has the dual functions of triglyceride hydrolase and ligand/bridging factor for receptor-mediated lipoprotein uptake.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.