



# Recombinant Glutamine synthetase 1 (glnA1)

<b>Product Code</b>	CSB-EP358639MVH
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0A591
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mycobacterium bovis (strain ATCC BAA-935 / AF2122/97)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTEKTPDDVF KLAKEKVEY VDVRFCDLPG IMQHFTIPAS AFDKSVFDDG LAFDGGSSIRG FQSIHESDML LLPDPETARI DPFRAAKTLN INFFVHDPFT LEPYSRDPRN IARKAENYLI STGIADTAYF GAEAEFYIFD SVSFDSTRANG SFYEVDASIG WWNTGAATEA DGSPNRGYKV RHKGGYFPVA PNDQYVDLRD KMLTNLINSG FILEKGGHHEV GSGGQAEINY QFNSSLHAAD DMQLYKYIIK NTAWQNGKTV TFMPKPLFGD NGSGMHCHQS LWKDGAPLMY DETYAGLSD TARHYIGLL HHAPSLLAFT NPTVNSYKRL VPGYEAPINL VYSQRNRSAC VRIPITG SNP KAKRLEFRSP DSSGNPYLAF SAMLMAGLDG IKNKIEPQAP VDKDLYELPP EEAASIPQTP TQLSDVIDRL EADHEYLTEG GVFTNDLIET WISFKRENEI EPVNIRPHPY EFALYYDV
<b>Source</b>	E.coli
<b>Target Names</b>	glnA1
<b>Protein Names</b>	Recommended name: Glutamine synthetase 1 EC= 6.3.1.2 Alternative name(s): Glutamate--ammonia ligase 1
<b>Expression Region</b>	1-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.