



# Recombinant Mouse Transcription factor MafF (Maff)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP013318MO-B                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Uniprot No.</b>       | O54791                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Immunogen Species</b> | Mus musculus (Mouse)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Sequence</b>          | MAVDPLSSKA LKVKRELSEN TPHLSDEALM GLSVRELNRN LRGLSAEEVT<br>RLKQRRRTLK NRGYAASCRV KRVCQKEELQ KQKSELEREV DKLARENAAM<br>RLELDALRGK CEALQGFARS VAAARGPAAL VAPASVITIV KSAPGPAPAA<br>DPAPCS                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Target Names</b>      | Maff                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Protein Names</b>     | Recommended name: Transcription factor MafF Alternative name(s): V-maf musculoaponeurotic fibrosarcoma oncogene homolog F                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Expression Region</b> | 1-156                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Protein Length</b>    | Full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Target Details</b>    | This protein is a basic leucine zipper (bZIP) transcription factor that lacks a transactivation domain. It is known to bind the US-2 DNA element in the promoter of the oxytocin receptor (OTR) gene and most likely heterodimerizes with other leucine zipper-containing proteins to enhance expression of the OTR gene during term pregnancy. The encoded protein can also form homodimers, and since it lacks a transactivation domain, the homodimer may act as a repressor of transcription. This gene may also be involved in the cellular stress response. Multiple transcript variants encoding two different isoforms have been found for this gene. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                                                                                   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                                                                                                                                                                                                |



of lyophilized form is 12 months at -20°C/-80°C.