



# Recombinant Amino-acid acetyltransferase (argA)

<b>Product Code</b>	CSB-EP358809EGX
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0A6C6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVKERKTELV EGFRHSVPYI NTHRGKTFVI MLGGEAIEHE NFSSIVNDIG LLHSLGIRLV VVYGARPQID ANLAAHHEP LYHKNI RVT D AKTLELVKQA AGTLQLDITA RLSMSLNNTPLQG AHINVVS GNFI AQPLG VDDGVDYCHS GRIRRIDEDA IHRQLDSGAI VLMGPVAVSV TGESFNLTSE EIATQLAIKL KAEKMIGFCS SQGV TND DGD IVSE LFPNEA QARVEAQEEK GDYNSGTVRF LRGAVKACRS GVRRCHLISY QEDGALLQEL FSRDGIGTQI VMESAEQIRR ATINDIGGIL ELIRPLEQQG ILVRRSREQL EMEIDKFTII QRDNTTIACA ALYPFPEEKI GEMACVAVHP DYRSSSRGEV LLERIAAQAK QSGLSKLFVL TTRSIHWFQE RGFTPVDIDL LPESKKQLYN YQRKSKVLMA DLG
<b>Source</b>	E.coli
<b>Target Names</b>	argA
<b>Protein Names</b>	Recommended name: Amino-acid acetyltransferase EC= 2.3.1.1 Alternative name(s): N-acetylglutamate synthase Short name= AGS Short name= NAGS
<b>Expression Region</b>	1-443
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.