



Recombinant Mouse Stromelysin-2 (Mmp10)

Product Code	CSB-EP014657MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O55123
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	F STFPGSPKWR KSHITYRIVN YTPDLPRQSV DSAIEKALKV WEEVTPLTFS RISEGEADIM ISFAVGEHGD FYPFDGPGQS LAHAYPPGPG FYGDVHFDDD EKWTLAPSGT NLFLVAAHEL GHSLGLFHSD KKESLMYPVY RFSTSPANFH LSQDDIEGIQ SLYGAGPSSD ATVPVLSVS PRPETPDKCD PALSFDSVST LRGEVLFFKD RYFWRRSHWN PEPEFHLISA FWPTLPDLD AAYEAHNTDS VLIFKGSQFW AVRGNEVQAG YPKGHTLGF PPTVKKIDAA VFEKEKKKTY FFVGDKYWRF DETRHVMDKG FPRQITDDFP GIEPQVDAVL HEFGFFYFFR GSSQFEFDPN ARTVTHILKS NSWLLC
Source	E.coli
Target Names	Mmp10
Protein Names	Recommended name: Stromelysin-2 Short name= SL-2 EC= 3.4.24.22 Alternative name(s): Matrix metalloproteinase-10 Short name= MMP-10 Transin-2
Expression Region	100-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP s are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. The enzyme encoded by this gene degrades proteoglycans and fibronectin. The gene is part of a cluster of MMP genes which localize to chromosome 11q22.3.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.