



# Recombinant Enolase (eno)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP358859EGX-B                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                  |
| <b>Uniprot No.</b>       | P0A6Q0                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Immunogen Species</b> | Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Sequence</b>          | SKIVKIIGR EIIDSRGNPT VEAEVHLEGG FVGMAAAPSG ASTGSREALE<br>LRDGDKSRFL GKGVTKAVAA VNGPIAQALI GKDAKDQAGI DKIMIDLDT<br>ENKSKFGANA ILAVSLANAK AAAAAKGMPL YEHIAELNGT PGKYSMPVPM<br>MNIINGGEHA DNNVDIQEFM IQPVGAKTVK EAIRMGSEVF HHLAKVLKAK<br>GMNTAVGDEG GYAPNLGSNA EALAVIAEAV KAAGYELGKD ITLAMDC AAS<br>EFYKDGKYVL AGEKNKAFTS EEFTHFLEEL TKQYPIVSIE DGLDESDWDG<br>FAYQTKVLGD KIQLVGDDLF VTNTKILKEG IEKGIANSIL IKFNQIGSLT<br>ETLAAIKMAK DAGYTAVISH RSGETEDATI ADLAVGTAAG QIKTGSMRSRS<br>DRVAKYNQLI RIEEALGEKA PYNGRKEIKG QA |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Target Names</b>      | eno                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Protein Names</b>     | Recommended name: Enolase EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 2-phosphoglycerate dehydratase                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Expression Region</b> | 2-432                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                         |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                  |