



# Recombinant Glycerol-3-phosphate dehydrogenase [NAD (P)+]

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP358874EOD   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | P0A6S9  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O157:H7  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MNQRNASMTV IGAGSYGTAL AITLARNGHE VVLWGHDPPEH IATLERDRCN<br>AAFLPDVPFP DTLHLESDLA TALAASRNIL VVPSHFVGE VLRQIKPLMR<br>PDARLVWATK GLEAETGRLL QDVAREALGD QIPLAVISGP TFAKELAAGL<br>PTAISLASTD QTFADDLQQL LHCCKSFRVY SNPDFIGVQL GGAVKNVIAI<br>GAGMSDGIGF GANARTALIT RGLAEMSRLG AALGADPATF MGMAGLGDLV<br>LTCTDNQSRN RRFGMMLGQG MDVQSAQEKI GQVVEGYRNT<br>KEVRELAHRF GVEMPITEEI YQVLYCGKNA REAALTLLGR ARKDERSSH      |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | gpsA  |
| <b>Protein Names</b>     | Recommended name: Glycerol-3-phosphate dehydrogenase [NAD(P)+] EC= 1.1.1.94 Alternative name(s): NAD(P)H-dependent glycerol-3-phosphate dehydrogenase   |
| <b>Expression Region</b> | 1-339   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |