



# Recombinant Glucose-1-phosphate adenylyltransferase (glgC)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP358880SZB   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | P0A6V4  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Shigella flexneri   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | VSLEKNDHL MLARQLPLKS VALILAGGRG TRKDLTNKR AKPAVHFGGK<br>FRIIDFALSN CINSGIRRMG VITQYQSHTL VQHIQRGWSF FNEEMNEFVD<br>LLPAQQRMKG ENWYRG TADA VTQNLDIIRR YKAEYVVILA GDHIYKQDYS<br>RMLIDHVEKG ARCTVACMPV PIEEASAFGV MAVDENDKII EFVEKPANPP<br>SMPNDPSKSL ASMGIVFDA DYLYELLEED DRDENS SHDF GKDLIPKITE<br>AGLAYAHPFP LSCVQSDPDA EPYWRDVGTL EAYWKANLDL ASVVPEDMY<br>DRNWPIRTYN ESLPPAKFVQ DRSGSHG MTL NSLVSGGCVI SGSVVVQSVL<br>FSRVRVNSFC NIDSAVLLPE VVVG RSCRLR RCVIDRACVI PEGMVIGENA<br>EEDARRFYRS EEGIVLV TRE MLRKLGHKQE R |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | glgC  |
| <b>Protein Names</b>     | Recommended name: Glucose-1-phosphate adenylyltransferase EC= 2.7.7.27<br>Alternative name(s): ADP-glucose pyrophosphorylase Short name= ADPGlc<br>PPase ADP-glucose synthase   |
| <b>Expression Region</b> | 2-431   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |