



Recombinant Adenylosuccinate synthetase (purA)

Product Code	CSB-EP358960EGX-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A7D5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	GNNVVVLGT QWGDEGKGI VDLLTERAKY VVRYQGGHNA GHTLVINGEK TVLHLIPSGI LRENVTSIIG NGVVLSPAAL MKEMKELEDR GIPVRERLLL SEACPLILDY HVALDNAREK ARGAKAIGTT GRGIGPAYED KVARRGLRVG DLFDKETFAE KLKEVMEYHN FQLVNYKAE AVDYQKVLDD TMAVADILTS MVVDVSDLLD QARQRGDFVM FEQAQGTLID IDHGTYPYVT SSNTTAGGVA TGSGLGPRYV DYVLGILKAY STRVGAGPFP TELFDETGEF LCKQGNEFGA TTGRRRRRTGW LDTVAVRRRAV QLNSLSGFCL TKLDVLDGLK EVKLCVAYRM PDGREVTTTP LAADDWKGVE PIYETMPGWS ESTFGVKDRS GLPQAALNYI KRIEELTGVP IDIISTGPDR TETMILRDPF DA
Source	E.coli
Target Names	purA
Protein Names	Recommended name: Adenylosuccinate synthetase Short name= AMPSase Short name= AdSS EC= 6.3.4.4 Alternative name(s): IMP--aspartate ligase
Expression Region	2-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.