



# Recombinant Human Sushi repeat-containing protein SRPX2 (SRPX2)

<b>Product Code</b>	CSB-YP022690HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O60687
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TWYAGSG YYPDESYNEV YAEVVPQAPA LDYRVPRWCY TLNIQDGEAT CYSKGGNYH SSLGTRCELS CDRGFRLIGR RSVQCLPSRR WSGTAYCRQM RCHALPFITS GTYTCTNGVL LDSRCDYSCS SGYHLEGDRS RICMEDGRWS GGEPVCVDID PPKIRCPHSR EKMAEPEKLT ARVYWDPLV KDSADGTITR VTLRGPEPGS HFPEGEHVIR YTAYDRAYNR ASCKFIVKVQ VRRCP TLKPP QHGYLTCTSA GDNYGATCEY HCDGGYDRQG TPSRVCQSSR QWSGSPICA PMKINVNVNS AAGLLDQFYE KQRLLIISAP DPSNRYYKMQ ISMLQQSTCG LDLRHVTIIE LVGQPPQEVG RIREQQLSAN IIEELRQFQR LTRSYFNMVL IDKQGIDRDR YMEPVTPEEI FTFIDYLLS NQELTQRREQ RDICE
<b>Source</b>	Yeast
<b>Target Names</b>	SRPX2
<b>Protein Names</b>	Recommended name: Sushi repeat-containing protein SRPX2
<b>Expression Region</b>	24-465
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.