



Recombinant Human PRAME family member 4 (PRAMEF4)

Product Code	CSB-MP018621HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O60810
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKMSIWTPPR LLELAGRSLL RDQALAMSTL EELPTELFPP LFMEAFSRRR CEALKLMVQS WPFRRPLRP LIKMPCLEAF QAVLDGLDAL LNLGVRPRRW KLQVLDLQDV CENFWMVWSE AMAHGCFLNA KRNKKPVEDC PRMRGRQPLT VFVELWLKNR TLDEYLTCLL LWVKQRKDLL HLCCCKLKIL GMPFRNIRSI LKMVNLDICIQ EVEVNCKWVL PILTQFTPYL GHMRNLQKLI LSHMDVSRVYV SPEQKKEIVT QFTTQFLKLR CLQKLYMNSV SFLEGHLDQL LSCLKTSLKF LTITNCVLE SDLKHL SQCP SISQLKTL DL SGIRLTNYSL VPLQILLEKV AATLEYLDLD DCGIIDSQVN AILPALSRFC ELNTFSFCGN PICMATLENL LSHTIILKNL CVELYPAPRE SYGADGTLCW SRFAQIRAEL MNRVRDLRHP KRILFCTDYC PDCGNRSFYD LEADQYCC
Source	Mammalian cell
Target Names	PRAMEF4
Protein Names	Recommended name: PRAME family member 4
Expression Region	1-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.