



Recombinant Escherichia coli Lipoprotein 34 (nlpB)

Product Code	CSB-MP359222ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A903
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	CSSDSR YKRQVSGDEA YLEAAPLAEL HAPAGMILPV TSGDYAIPVT NGSGAVGKAL DIRPPAQPLA LVSGARTQFT GDTASLLVEN GRGNTLWPQV VSVLQAKNYT ITQRDDAGQT LTTDWWQWNR LDEDEQYRGR YQISVKPQGY QQAVTVKLLN LEQAGKPVAD AASMQRYS TE MMNVISAGLD KSATDAANAA QNRASTTMDV QSAADDTGLP MLVVRGPFNV VWQRLPAALE KVG MKVTDST RSQGNMAVTY KPLSDSDWQE LGASDPGLAS GDYKLQVGDL DNRSSLQFID PKGHTLTQSQ NDALVAVFQA AFSK
Source	Mammalian cell
Target Names	bamC
Protein Names	Recommended name: Lipoprotein 34
Expression Region	25-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.