



# Recombinant Anaerobic glycerol-3-phosphate dehydrogenase subunit A (glpA)

<b>Product Code</b>	CSB-EP359277EOD-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0A9C1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKTRDSQSSD VIIIGGGATG AGIARDCALR GLRVILVERH DIATGATGRN HGLLHSGARY AVTDAESARE CISENQILKR IARHCVEPTN GLFITLPEDD LSFQATFIRA CEEAGISAEA IDPQQARIIE PAVNPALIGA VKVPDGTVDP FRLTAANMLD AKEHGAVILT AHEVTGLIRE GATVCGVRVR NHLTGETQAL HAPVVVNAAG IWGQHIAEYA DLRIRMFPK GSELLIMDHRI NQHVINRCRK PSDADILVPG DTISLIGTTS LRIDYNEIDD NRVTAAEVDI LLREGEKLAP VMAKTRILRA YSGVRPLVAS DDDPSGRNVS RGIVLLDHAERDGLDGFITI TGGKLMYRL MAEWATDAVC RKLGNTRPCT TADLALPGSQ EPAEVTLRKV ISLPAPLRGS AVYRHGDRTP AWLSEGRLHR SLVCECEAVT AGEVQYAVEN LNVNSLLDLR RRTRVGMGTC QGELCACRAA GLLQRFNVTT SAQSIEQLST FLNERWKGVQ PIAWGDALRE SEFTRWVYQG LCGLEKEQKD AL
<b>Source</b>	E.coli
<b>Target Names</b>	glpA
<b>Protein Names</b>	Recommended name: Anaerobic glycerol-3-phosphate dehydrogenase subunit A Short name= G-3-P dehydrogenase EC= 1.1.5.3
<b>Expression Region</b>	1-542
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.