



Recombinant Escherichia coli HTH-type transcriptional regulator metR (metR)

Product Code	CSB-BP359291ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A9F9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MIEVKHLKTL QALRNCGLA AAAATLHQTQ SALSHQFSDL EQRLGFRLFV RKSQPLRFTP QGEILLQLAN QVLPQISQAL QACNEPQQTR LRIAIECHSC IQWLTPALEN FHKNWPQVEM DFKSGVTFDP QPALQQGELD LVMTSDILPR SGLHYSPMFD YEVRVLVLPD HPLAAKTRIT PEDLASELL IYPVQRSRLD VWRHFLQPAG VSPSLKSVDN TLLLIQMVA A RMGIAALPHW VVESFERQGL VVTKTLGEGE WSRLYAAVRD GEQRQPVTEA FIRSARNHAC DHLPFVKSAE RPTYDAPTVR PGSPARL
Source	Baculovirus
Target Names	metR
Protein Names	Recommended name: HTH-type transcriptional regulator metR
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.