



Recombinant Cyclopropane-fatty-acyl-phospholipid synthase (cfa)

Product Code	CSB-YP359300EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A9H8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	SSSCIEEVS VPDDNWYRIA NELLSRAGIA INGSAPADIR VKNPDDFFKRV LQEGSLGLGE SYMDGWWECD RLDMFFSKVL RAGLENQLPH HFKDTLRIAG ARLFNLQSKK RAWIVGKEHY DLGNDLFSRM LDPFMQYSCA YWKDADNLES AQQAKLMIC EKLQLKPGMR VLDIGCGWGG LAHYMASNYD VSVVGVITISA EQQKMAQERC EGLDVTILLQ DYRDLNDQFD RIVSVGMFEH VGPKNYDTYF AVVDRNLKPE GIFLLHTIGS KKTDLNVDPW INKYIFPNGC LPSVRQIAQS SEPHFVMEDW HNFADYDTT LMAWYERFLA AWPEIADNYS ERFKRMFTYY LNACAGAFRA RDIQLWQVVF SRGVENGLRV AR
Source	Yeast
Target Names	cfa
Protein Names	Recommended name: Cyclopropane-fatty-acyl-phospholipid synthase Short name= CFA synthase Short name= Cyclopropane fatty acid synthase EC= 2.1.1.79
Expression Region	2-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.