



# Recombinant Ribonuclease G (rng)

<b>Product Code</b>	CSB-MP359307EGX
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0A9J1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TAELLVNVT PSETRVAYID GGILQEIHIE REARRGIVGN IYKGRVSRVL PGMQAAFVDI GLDKAAFLHA SDIMPHTECV AGEEQKQFTV RDISELVRQG QDLMVQVVKD PLGTKGARLT TDITLPSRYL VFMPGASHVG VSQRIESESE RERLKKVVAE YCDEQGGFII RTAAEGVGEA ELASDAAYLK RVWTKVMERK KRPQTRYQLY GELALAQRVL RDFADAELDR IRVDSRLTYE ALLEFTSEYI PEMTSKLEHY TGRQPIFDLF DVENEIQRAL ERKVELKSGG YLIIDQTEAM TTVDINTGAF VGHRNLDDTI FNTNIEATQA IARQLRLRNL GGIIIDFID MNNEHDRRRV LHSLEQALS KDRVKTSVNGF SALGLVEMTR KRTRESIEHV LCNECPTCHG RGTVKTIVETV CYEIMREIVR VHHAYDSDRF LUYASPAVAE ALKGEESHSL AVEVEIFVGKQ VKVQIEPLYN QEQFDVMM
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rng
<b>Protein Names</b>	Recommended name: Ribonuclease G Short name= RNase G EC= 3.1.26.- Alternative name(s): Cytoplasmic axial filament protein
<b>Expression Region</b>	2-489
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.