



Recombinant Escherichia coli PhoH-like protein (ybeZ)

Product Code	CSB-BP359312ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A9K3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MNIDTREITL EPADNARLLS LCGPFDDNIK QLERRLGIEI NRRDNHFKLT GRPICVTAAA DILRSLYVDT APMRGQIQDI EPEQIHLAIK EARVLEQSAE SVPEY GKAVN IKTKR GVIKP RTPNQAQYIA NILDHDITFG VGPAGTGKTY LAVAAVDAL ERQEIRRILL TRPAVEAGEK LGFLPGDLSQ KVDPYLRPLY DALFEMLGFE KVEKLIERNV IEVAPLAYMR GRTLNDAFII LDESQNTTIE QMKMFLTRIG FNSKAVITGD VTQIDLPRNT KSGLRHAIEV LADVEEISFN FFHSE DVVRH PVVARIVNAY EAWEEAEQKR KAALAAERKR EEQE QK
Source	Baculovirus
Target Names	ybeZ
Protein Names	Recommended name: PhoH-like protein
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.