



Recombinant Uncharacterized ABC transporter ATP-binding protein YbbA (ybbA)

Product Code	CSB-BP359345SZB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A9U0
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	≥85% (SDS-PAGE)
Sequence	MPAENIVEVH HLKKSVMGQGE HELSILTGVE LVVKRGETIA LVGESGSGKS TLLAILAGLD DGSSGEVSLV GQPLHNMDEE ARAKLRAKHV GFVFQSFMLI PTLNALLENVE LPALLRGESS AESRNGAKAL LEQLGLGKRL DHLPAQLSGG EQQRVALARA FNGRPDVLFA DEPTGNLDRQ TGDKIADLLF SLNREHGTTL IMVTHDLQLA ARCDRCLRLV NGQLQEEA
Source	Baculovirus
Target Names	ybbA
Protein Names	Recommended name: Uncharacterized ABC transporter ATP-binding protein YbbA
Expression Region	1-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.