



Recombinant Uncharacterized ABC transporter ATP-binding protein YbiT (ybiT)

Product Code	CSB-BP359347EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A9U4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MLVSSNVTMQ FGSKPLFENI SVKFGGGNRY GLIGANGSGK STFMKILGGD LEPTLGNVSL DPNERIGKLR QDQFAFEFT VLDTVIMGHK ELWEVKQERD RIYALPEMSE EDGYKVADLE VKYGEMDGYS AEARAGELLL GVGIPVEQHY GPMSEVAPGW KLRVLLAQAL FADPDILLDD EPTNNLDIDT IRWLEQVLNE RDSTMIISH DRHFLNMVCT HMAADLDY GEL RVYPGNYDEY MTAATQARER LLADNAKKKA QIAELQSFVS RFSANASKSR QATSRARQID KIKLEEVKAS SRQNPFI RFE QDKKLF RNAL EVEGLTKGFD NGPLFKNLNL LLEVGEK LAV LGTNGVGKST LLKTLV GDLQ PDSGTVKWSE NARIGYYAQD HEYEFENDLT VFEWMSQWKQ EGDDEQAVRS ILGRLLFSQD DIKKPAKVL S GGEKGRMLFG KLMMQKPNIL IMDEPTNHLD MESIESLNMA LELYQGTLIF VSHDREFVSS LATRILEITP ERVIDFSGNY EDYLRSKGIE
Source	Baculovirus
Target Names	ybiT
Protein Names	Recommended name: Uncharacterized ABC transporter ATP-binding protein YbiT
Expression Region	1-530
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.