



Recombinant Uncharacterized lipoprotein yeaY (yeaY)

Product Code	CSB-BP359420EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AA92
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	CVTVPDAI KGSSPTPQQD LVRVMSAPQL YVGQEARFGG KVVAVQNQQG KTRLEIATVP LDSGARPTLG EPSRGRIYAD VNGFLDPVDF RGQLVTVVGP ITGAVDGKIG NTPYKFMVMQ VTGYKRWHLT QQVIMPPQPI DPWFYGGRGW PYGYGGWGWY NPGPARVQTV VTE
Source	Baculovirus
Target Names	yeaY
Protein Names	Recommended name: Uncharacterized lipoprotein yeaY
Expression Region	23-193
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.