



Recombinant 3-octaprenyl-4-hydroxybenzoate carboxy-lyase (ubiD)

Product Code	CSB-BP359432EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AAB5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MDAMKYNDLR DFLTLEQQG ELKRITLPVD PHLEITEIAD RTLRAGGPAL LFENPKGYSM PVLCLNLFQTP KRVAMGMGQE DVSALREVGK LLAFLKEPEP PKGFRDLFDK LPQFKQVLNM PTKRLRGAPC QQKIVSGDDV DLNRIPIMTC WPEDAAPLIT WGLTVTRGPH KERQNLGIYR QQLIGKNKLI MRWLSHRGGGA LDYQEWCAAH PGERFPVSVVA LGADPATILG AVTPVPDTLS EYAFAGLLRG TKTEVVKCIS NDLEVPASAE IVLEGYIEQG ETAPEGPYGD HTGYINEVDS FPVFTVTHIT QREDAIYHST YTGRPPDEPA VLGVALNEVF VPILQKQFPE IVDFYLPPEG CSYRLAVVTI KKQYAGHAKR VMMGVWSFLR QFMYTKFVIV CDDVDNARDW NDVIWAITTR MDPARDTVLV ENTPIDYLDF ASPVSGLGSK MGLDATNKWP GETQREWGRP IKKDPDVVAH IDAIWDELAI FNNGKSA
Source	Baculovirus
Target Names	ubiD
Protein Names	Recommended name: 3-octaprenyl-4-hydroxybenzoate carboxy-lyase EC=4.1.1.- Alternative name(s): Polyprenyl p-hydroxybenzoate decarboxylase
Expression Region	1-497
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.