



Recombinant Universal stress protein E (uspE)

Product Code	CSB-EP359434EGX-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AAC1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	AMYQNMLVV IDPNQDDQPA LRRAVYLHQR IGGKIKAFLLP IYDFSYEMTT LLSPDERTAM RQGVISQRTA WIHEQAKYYL NAGVPIEIKV VWHNRPFEAI IQEVISGGHD LVLKMAHQHD RLEAVIFTPT DWHLLRKCPSPVWVMVKDQPW PEGGKALVAV NLASEEPYHN ALNEKLVKET IELAEQVNHT EVHLVGAYPV TPINIAIELP EFDPSVYNDA IRGQHLLAMK ALRQKFGINE NMTHVEKGLP EEVIPDLAEH LQAGIVVLGT VGRTGISAAF LGNTAEQVID HLRCDLLVIK PDQYQTPVEL DDEEDD
Source	E.coli
Target Names	uspE
Protein Names	Recommended name: Universal stress protein E
Expression Region	2-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.