



Recombinant Escherichia coli Probable anaerobic dimethyl sulfoxide reductase chain ynfG (ynfG)

Product Code	CSB-BP359469ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AAJ1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	TTQYGFFID SSRCTGCKTC ELACKDFKDL GPEVSFRRYIY EYAGGDWQED NGVWHQNVFA YYLSISCNHC DDPACTKVCP SGAMHKREDG FVVVDEDVCI GCRYCHMACP YGAPQYNAEK GHMTKCDGCV SRVAEGKQPI CVESCPLRAL EFGPIELRQ KHGTLAAVAP LPRAHFTKPN IVIKPNANSR PTGDTTGYLA NPEEV
Source	Baculovirus
Target Names	ynfG
Protein Names	Recommended name: Probable anaerobic dimethyl sulfoxide reductase chain ynfG Alternative name(s): DMSO reductase iron-sulfur subunit ynfG
Expression Region	2-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.