



Recombinant Escherichia coli Replication-associated recombination protein A (rarA)

Product Code	CSB-YP359545ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AAZ4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSNLSLDFSD NTFQPLAARM RPENLAQYIG QQHLLAAGKP LPRAIEAGHL HSMILWGPPG TGKTTLAEVI ARYANADVER ISAVTSGVKE IREAIERARQ NRNAGRRTIL FVDEVHRFNK SQQDAFLPHI EDGTITFIGA TTENPSFELN SALLSRARVY LLKSLSTEDI EQVLTQAMED KTRGYGGQDI VLPDETRRAI AELVNGDARR ALNTLEMMAD MAEVDSDGKR VLKPELLTEI AGERSARFDN KGDRFYDLIS ALHKSVRGS A PDAALYWYAR IITAGGDPLY VARRCLAIAS EDVGNADPRA MQVAIAAWDC FTRVGPAEGE RAIQAIVYL ACAPKSNAVY TAFKAALADA RERPDYDVPV HLRNAPTKLM KEMGYGQEYR YAHDEANAYA AGEVYFPPEI AQTRYFPTN RGLEGKIGEK LAWLAEQDQN SPIKRYR
Source	Yeast
Target Names	rarA
Protein Names	Recommended name: Replication-associated recombination protein A
Expression Region	1-447
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.