



Recombinant Putative 8-amino-7-oxononanoate synthase/2-amino-3-ketobutyrate coenzyme A ligase (Z5044, ECs4495)

Product Code	CSB-BP359586EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AB78
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MRGEFYQQLT NDLETARAEG LFKEERIITS AQQADITVAD GSHVINFCAN NYLGLANHPD LIAAAKAGMD SHGFGMASVR FICGTQDSHK ELEQKLA AFL GMEDAILYSS CFDANGGLFE TLLGAEDAI SDALNHASII DGVRLCKAKR YRYANNDMQE LEARLKEARE AGARHVLIAT DGVFSMDGVI ANLKGVCDLA DKYDALVMVD DSHAVGFVGE NGRGSHEYCD VMGRVDIITG TLGKALGGAS GGYTAAARKEV VEWLRQRSRP YLFSNSLAPA IVAASIKVLE MVEAGSELRD RLWANARQFR EQMSAAGFTL AGADHAIIPV MLGDAVVAQK FARELQKEGI YVTGFFYPVV PKGQARIRTQ MSAAH TPEQI TRAVEAFTRI GKQLGVIA
Source	Baculovirus
Target Names	kbl
Protein Names	Recommended name: Putative 8-amino-7-oxononanoate synthase/2-amino-3-ketobutyrate coenzyme A ligase Short name= AONS/AKB ligase EC= 2.3.1.29 EC= 2.3.1.47 Alternative name(s): 7-keto-8-amino-pelargonic acid synthase Short n
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.