



Recombinant Thiamine biosynthesis lipoprotein ApbE (apbE)

Product Code	CSB-EP359588EOD-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AB86
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	C DQKPQPAKTH ATEVTVLEGK TMGTFWRASI PGIDAKRSAE LKEKIQTQLD ADDQLLSTYK KDSALMRFND SQSLSPWPVS EAMADIVTTS LRIGAKTDGA MDITVGPLVN LWGFGPEQQP VQIPSQEQID AMKAKTGLQH LTVINQSHQQ YLQKDLPDLY VDLSTVGEY AADHLARLME QEGISRYLVS VGGALNSRGM NGEGLPWRVA IQKPTDKENA VQAVVDINGH GISTSGSYRN YYELDGKRLS HVIDPQTGRP IEHNLVSVTV IAPTALEADA WDTGLMVLGP EKAKEVVRRE GLAVYMITKE GDSFKTWMSP QFKSFLVSEK N
Source	E.coli
Target Names	apbE
Protein Names	Recommended name: Thiamine biosynthesis lipoprotein ApbE
Expression Region	20-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.