



Recombinant Human WD repeat-containing protein 1 (WDR1)

Product Code	CSB-BP025995HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O75083
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPYEIKKVFASLPQVERGVSKIIGGDPKGNNFLYTNGKCVILRNIDNPALADIYT EHAHQVVVAKYAPSGFYIASGDVSGKLRWDTTQKEHLLKYEYQPFAGKIKDIA WTEDSKRIAVVGEGREKFGAVFLWDSGSSVGEITGHNVINSVDIKQSRPYRL ATGSDDNCAAFFEGPPFKFKFTIGDHSRFVNCVRFSPDGNRFATASADGQIYI YDGKTGEKVCALGGSKAHDGGIYAIWSPDSTHLLSASGDKTSKIWDVSVNSV VSTFPMGSTVLDQQLGCLWQKDHLSSVLSGYINYLDRNNPSKPLHVIKGHSK SIQCLTVHKNGGKSYIYSGSHDGHINYWDESETGENDSFAGKGHTNQVSRMTV DESGQLISCSMDDTVRYTSLMLRDYSGQGQVVKLDVQPKCVAVGPGGYAVVV CIGQIVLLKDQRKCFSIDNPGYEPEVVAVHPGGDTVAIGGVDGNVRLYSILGTT LKDEGKLLKAKGPVTDVAYSHDGAFLAVCDASKVVTVFSVADGYSENNVFIG HHAKIVCLAWSPDNEHFASGGMDMMVYVWTLSDPETRVKIQDAHRLHHVSSL AWLDEHTLVTTSHDASVKEWTITY
Source	Baculovirus
Target Names	WDR1
Protein Names	Recommended name: WD repeat-containing protein 1 Alternative name(s): Actin-interacting protein 1 Short name= AIP1 NORI-1
Expression Region	1-606aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.