



Recombinant Human Zinc finger protein 253 (ZNF253)

Product Code	CSB-BP026613HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75346
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGPLQFRDVA IEFSLLEWHC LDTAQRNLYR DVMLENYRNL VFLGIVVSKP DLVTCLEQGK KPLTMERHEM IAKPPVMSSH FAQDLWPENI QNSFQIGMLR RYEECRHDNL QLKKGCKSVG EHKVHKGGYN GLNQCLTTTQ KEIFQCCKYK KVFHKFSNSN TYKTRHTGIN LFKCIICGKA FKRSSTLTTH KKIHTGEKPY RCEECGKAFN QSANLTTHKR IHTGEKPYRC EECGKAFKQS SNLTTHKKIH TGEKPYKCEE CGKAFNRSTD LTTHKIVHTG EKPYKCEECEG KAFKHPSHVT THKKIHTRGK PYNCEECEGKS FKHCNSLTIH KRIHTGEKPY KCEECGKAFH LSSHLTTHKI LHTGEKPYRC RECGKAFNHS TTLFSHEKIH TGEKPYKCDE CGKTFTWPSI LSKHKRHTG EKPYKCEECEG KSFTASSTLT THKRIHTGEK PYKCEECEGKA FNWSSDLNKH KKIHIERKPY IVKNVTDLLN VPPLLSIR
Source	Baculovirus
Target Names	ZNF253
Protein Names	Recommended name: Zinc finger protein 253 Alternative name(s): Bone marrow zinc finger 1 Short name= BMZF-1 Zinc finger protein 411
Expression Region	1-499
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.