



# Recombinant Human Zinc finger protein 253 (ZNF253)

<b>Product Code</b>	CSB-MP026613HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O75346
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGPLQFRDVA IEFSLLEWHC LDTAQRNLYR DVMLENYRNL VFLGIVVSKP DLVTCLEQGK KPLTMERHEM IAKPPVMSSH FAQDLWPENI QNSFQIGMLR RYEECRHDNL QLKKGCKSVG EHKVHKGGYN GLNQCLTTTQ KEIFQCCKYG KVFHKFSNSN TYKTRHTGIN LFKCIICGKA FKRSSTLTTH KKIHTGEKPY RCEECGKAFN QSANLTTHKR IHTGEKPYRC EECGKAFKQS SNLTTHKKIH TGEKPYKCEE CGKAFNRSTD LTTHKIVHTG EKPYKCEECG KAFKHPSHVT THKKIHTRGK PYNCEECGKS FKHCNSLTIH KRIHTGEKPY KCEECGKAFH LSSHLTTHKI LHTGEKPYRC RECGKAFNHS TTLFSHEKIH TGEKPYKCDE CGKTFTWPSI LSKHKRHTG EKPYKCEECG KSFTASSTLT THKRIHTGEK PYKCEECGKA FNWSSDLNKH KKIHIERKPY IVKNVTDLLN VPPLLSIR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF253
<b>Protein Names</b>	Recommended name: Zinc finger protein 253 Alternative name(s): Bone marrow zinc finger 1 Short name= BMZF-1 Zinc finger protein 411
<b>Expression Region</b>	1-499
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.