



Recombinant Protein mpaA (mpaA)

Product Code	CSB-EP359839EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0ACV7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MTVTRPRAER GAFPPGTEHY GRSLLGAPLI WFPAPAASRE SGLILAGTHG DENSSVVTLS CALRTLTPSL RRHHVVLNVN PDGCQLGLRA NANGVDLNRN FPAANWKEGE TVYRWNSAAE ERDVLLTGD KPGSEPETQA LCQLIHRIQP AWVVSFHDPL ACIEDPRHSE LGEWLAQAFE LPLVTSVGYE TPGSFGSWCA DLNLHCITAE FPPISSDEAS EKYL FAMANL LRWHPKDAIR PS
Source	E.coli
Target Names	mpaA
Protein Names	Recommended name: Protein mpaA
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.