



# Recombinant Uncharacterized protein yeaG (yeaG)

<b>Product Code</b>	CSB-YP359853EOD
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0ACY5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNIFDHYRQR YEAAKDEEFT LQEFLTTCRQ DRSAYANAAE RLLMAIGEPV MVDTAQEPRL SRLFSNRVIA RYPAFEFYG MEDAIEQIVS YLKHAAQGLE EKKQILYLLG PVGGGKSSLA ERLKSLMQLV PIYVLSANGE RSPVNDHPFC LFNPQEDAQI LEKEYGIPRR YLGTIMSPWA AKRLHEFGGD ITKFRVVKVW PSILQQIAIA KTEPGDENNQ DISALVGKVD IRKLEHYAQN DPDAYGYSGA LCRANQGIME FVEMFKAPIK VLHPLLTATQ EGNYNGTEGI SALPFNGIIL AHSNESEWVT FRNNKNNEAF LDRVYIVKVP YCLRISEEIK IYEKLLNHSE LTHAPCAPGT LETLSRFSIL SRLKEPENSS IYSKMRVYDG ESLKDTDPKA KSYQEYRDYA GVDEGMNGLS TRFAFKILSR VFNFDHVEVA ANPVHLFYVL EQQIEREQFP QEQAERYLEF LKGYLIPKYA EFIGKEIQTA YLESYSEYGG NIFDRYVTYA DFWIQDQEYR DPDTGQLFDR ESLNAELEKI EKPAGISNPK DFRNEIVNFV LRARANNSGR NPNWTSYEKL RTVIEKKMFS NTEELLPVIS FNAKTSTDEQ KKHDDFVDRM MEKGYTRKQV RLLCEWYLRV RKSS
<b>Source</b>	Yeast
<b>Target Names</b>	yeaG
<b>Protein Names</b>	Recommended name: Uncharacterized protein yeaG
<b>Expression Region</b>	1-644
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.