



Recombinant Human Deformed epidermal autoregulatory factor 1 homolog (DEAF1)

Product Code	CSB-EP006644HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75398
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEDSDSAAKQ LGLAEAAVA AAAAVAAAAA AAAGGEAEEP VLSRDEDSEE DADSEAERET PRVTAVAVMA AEPGHMDMGA EALPGPDEAA AAAFAEVTT VTVANVGAAA DNVFTTSVAN AASISGHVLS GRTALQIGDS LNTEKATLIV VHTDGSIVET TGLKGPAAPL TPGPQSPPTP LAPGQEKGGT KYNWDPSVYD SELPVRCRNI SGTLYKNRLG SGGRGRCIKQ GENWYSPTEF EAMAGRASSK DWKRSIRYAG RPLQCLIQDG ILNPHAASCT CAACCCDDMTL SGPVRLFVPY KRRKKENELP TTPVKKDSPK NITLLPATAA TTFTVTPSGQ ITTSGALTFD RASTVEATAV ISESPAQGDV FAGATVQEAS VQPPCRASHP EPHYPGYQDS CQIAPFPEAA LPTSHPKIVL TSLPALAVPP PTPTKAAPPA LVNGLELSEP RSWLYLEEMV NSLLNTAQQK KTLFEQAKHA STYREAAATNQ AKIHADAERK EQSCVNCGRE AMSECTGCHK VNYCSTFCQR KDWKDHQHIC GQSAAVTVQA DEVHVAESVM EKVTV
Source	E.coli
Target Names	DEAF1
Protein Names	Recommended name: Deformed epidermal autoregulatory factor 1 homolog Alternative name(s): Nuclear DEAF-1-related transcriptional regulator Short name= NUDR Suppressin Zinc finger MYND domain-containing protein 5
Expression Region	1-565
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.