



# Recombinant Human Deformed epidermal autoregulatory factor 1 homolog (DEAF1)

<b>Product Code</b>	CSB-YP006644HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O75398
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEDSDSAAKQ LGLAEAAVA AAAAVAAAAA AAAGGEAEEP VLSRDEDSEE DADSEAERET PRVTAVAVMA AEPGHMDMGA EALPGPDEAA AAAFAEVTT VTVANVGAAA DNVFTTSVAN AASISGHVLS GRTALQIGDS LNTEKATLIV VHTDGSIVET TGLKGPAAPL TPGPQSPPTP LAPGQEKGGT KYNWDPSVYD SELPVRCRNI SGTLYKNRLG SGGRGRCIKQ GENWYSPTEF EAMAGRASSK DWKRSIRYAG RPLQCLIQDG ILNPHAASCT CAACCCDDMTL SGPVRLFVPY KRRKKENELP TTPVKKDSPK NITLLPATAA TTFTVTPSGQ ITTSGALTFD RASTVEATAV ISESPAQGDV FAGATVQEAS VQPPCRASHP EPHYPGYQDS CQIAPFPEAA LPTSHPKIVL TSLPALAVPP PTPTKAAPPA LVNGLELSEP RSWLYLEEMV NSLLNTAQQK KTLFEQAKHA STYREAATNQ AKIHADAERK EQSCVNCGRE AMSECTGCHK VNYCSTFCQR KDWKDHQHIC GQSAAVTVQA DEVHVAESVM EKVTV
<b>Source</b>	Yeast
<b>Target Names</b>	DEAF1
<b>Protein Names</b>	Recommended name: Deformed epidermal autoregulatory factor 1 homolog Alternative name(s): Nuclear DEAF-1-related transcriptional regulator Short name= NUDR Suppressin Zinc finger MYND domain-containing protein 5
<b>Expression Region</b>	1-565
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.