



Recombinant D-alanyl-D-alanine-carboxypeptidase/endopeptidase AmpH (ampH)

Product Code	CSB-EP359897EOD-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AD71
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	AQPITEPEF ASDIVDRYAD HIFYGSGATG MALVVIDGNQ RVFRSYGETR PGNNVRPQLD SVVRIASLTK LMTSEMLVKL LDQGTVKLND PLSKYAPPGA RVPTYNGTPI TLVNLATHTS ALPREQPGGA AHRPVFVWPT REQRWKYLST AKLKAAPGSQ AAYSNLAFDL LADALANASG KPYTQLFEEQ ITRPLGMKDT TYTPSPDQCR RLMVAERGAS PCNNTLAAIG SGGVYSTPGD MMRWMQQYLS SDFYQRSNQA DRMQTLIYQR AQFTKVI GMD VPGKADALGL GWVYMAPKEG RPGIIQKTGG GGGFITYMAM IPQKNIGAFV VVTRSPLTRF KNMSDGINDL VTELSGNKPL VIPAS
Source	E.coli
Target Names	ampH
Protein Names	Recommended name: D-alanyl-D-alanine-carboxypeptidase/endopeptidase AmpH EC= 3.4.-.- Alternative name(s): DD-alanine-endopeptidase DD-carboxypeptidase Penicillin-binding protein Short name= PBP
Expression Region	22-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.