



Recombinant Human Transcription factor Maf (MAF)

Product Code	CSB-EP013312HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75444
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MASELAMSNS DLPTSPLAME YVNDFDLMKF EVKKEPVETD RIISQCGRLI AGGSLSTPM STPCSSVPPS PSFSAPSPGS GSEQKAHLED YYWMTGYPPQ LNPEALGFSP EDAVEALISN SHQLQGGFDG YARGAQQALAA AAGAGAGASL GGSGEEMGPA AAVVSAVIAA AAAQSGAGPH YHHHHHHAAG HHHHPTAGAP GAAGSAAASA GGAGGAGGGG PASAGGGGGG GGGGGGGGAA GAGGALPHHH AAGGLHFDDR FSDEQLVTMS VRELNRQLRG VSKEEVIRLK QKRRTLKNRG YAQSCRFKRV QQRHVLESEK NQLLQQVDHL KQEISRLVRE RDAYKEYEK LVSSGFRENG SSSDNPSSPE FFM</p>
Source	E.coli
Target Names	MAF
Protein Names	Recommended name: Transcription factor Maf Alternative name(s): Proto-oncogene c-Maf V-maf musculoaponeurotic fibrosarcoma oncogene homolog
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.