



Recombinant Human Zinc finger protein 324A (ZNF324)

Product Code	CSB-EP026673HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75467
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAFEDVAVYF SQEEWGLLDT AQRALYRRVM LDNFALVASL GLSTSRPRVV IQLERGEEPW VPSGTDTTLS RTTYRRRNPQ SWSLTEDRDV SGEWPRAFPD TPPGMTTSVF PVAGACHSVK SLQRQRGASP SRERKPTGVS VIYWERLLLG SGSGQASVSL RLTSPLRPE GVRLREKTLT EHALLGRQPR TPERQKPCAQ EVPGRFTGSA QDLEAAGGRG HHRMGAVWQE PHRLGGQEP STWDELGEAL HAGEKSFECR ACSVFVKSS DLLKHLRHT GERPYECAQC GKAFSQTSHL TQHQRHSGE TPYACPVCVK AFRHSSSLVR HQRIHTAEKS FRCSECGKAF SHGSNLSQHR KIHAGGRPYA CAQCGRRFCR NSHLIQHERT HTGEKPFVCA LCGAAFSQGS SLFKHQRVHT GEKPFACPQC GRAFSHSSNL TQHQLLHTGE RPFRCVDCGK AFAKGAVLLS HRIHTGEKP FVCTQCGRF RERPALFHHQ RIHTGEKTVR RSRASLHPQA RSVAGASSEG APAKETEPTP ASGPAAVSQP AEV
Source	E.coli
Target Names	ZNF324
Protein Names	Recommended name: Zinc finger protein 324A Alternative name(s): Zinc finger protein ZF5128
Expression Region	1-553
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.