



Recombinant Human Transcription initiation protein SPT3 homolog (SUPT3H)

Product Code	CSB-MP022956HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75486
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<pre>MYSRGRSQGR GTAEATANSP SPPIAPSHSR VTFSLSTLHT LSPPPRPFPS VSRAAAQKPH HLHPHILLAG SAAVPPRVLK AEMNNTAASP MSTATSSSSGR STGKSISFAT ELQSMMSYLG DARRPLHETA VLVEDVVHTQ LINLLQQAEE VSQLRGARVI TPEDLLFLMR KDKKKLRLL KYMFIRDYKS KIVKGIDEDD LLEDKLSGSN NANKRQKIAQ DFLNSIDQTG ELLAMFEDDE IDEVKQERME RAERQTRIMD SAQYAEFCES RQLSFSKKAS KFRDWLDCSS MEIKPNVVAM EILAYLAYET VAQLVDLALL VRQDMVTKAG DPFSHAISAT FIQYHNSAES TAACGVEAHS DAIQPCHIRE AIRRYSHRIG PLSPFTNAYR RNGMAFLAC</pre>
Source	Mammalian cell
Target Names	SUPT3H
Protein Names	Recommended name: Transcription initiation protein SPT3 homolog Alternative name(s): SPT3-like protein
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.