



Recombinant Lipopolysaccharide export system protein IptA (IptA)

Product Code	CSB-MP360011EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0ADV2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	VTG DTDQPIHIES DQQSLDMQGN VVTFTGNVIV TQGTIKINAD KVVVTRPGGE QGKEVIDGYG KPATFYQMVD NGKPVEGHAS QMHYELAKDF VVLTGNAYLQ QVDSNIKGDK ITYLVKEQKM QAFSDKGKRV TTVLVPSQLQ DKNNKGQTPA QKKGN
Source	Mammalian cell
Target Names	IptA
Protein Names	Recommended name: Lipopolysaccharide export system protein IptA
Expression Region	28-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.