



Recombinant Uncharacterized protein yhbW (yhbW)

Product Code	CSB-MP360013EGX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0ADV6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MTDKTIAFSL LDLAPIPEGS SAREAFSHSL DLARLAEKRG YHRYWLAEEH NMTGIASAAT SVLIGYLAAN TTTLHLGSGG VMLPNHSPLV IAEQFGTLNT LYPGRIDLGL GRAPGSDQRT MMALRRHMSG DIDNFPRDVA ELVDWFDARD PNPHVRPVPY YGEKIPVWLL GSSLYSAQLA AQLGLPFAPA SHFAPDMLFQ ALHLYRSNFK PSARLEKPYA MVCINIIAAD SNRDAEFLFT SMQQAFVKLR RGETGQLPPP IQNMDQFWSP SEQYGVQQAL SMSLVGDKAK VRHGLQSILR ETDADEIMVN GQIFDHQARL HSFELAMDVK EELLG
Source	Mammalian cell
Target Names	yhbW
Protein Names	Recommended name: Uncharacterized protein yhbW
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.