



Recombinant Escherichia coli Multidrug efflux pump subunit AcrA (acrA)

Product Code	CSB-YP360036ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0AE06
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	CDDKQA QQGGQQMPAV GVVTVKTEPL QITTELPGR T SAYRIAEVRP QVSGIILKRN FKEGSDIEAG VSLYQIDPAT YQATYDSAKG DLAKAQAAAN IAQLTVNRYQ KLLGTQYISK QEYDQALADA QQANA AVTAA KAAVETARIN LAYTKVTSPI SGRIGKSNVT EGALVQNGQA TALATVQQLD PIYVDVTQSS NDFLRLKQEL ANGLTKQENG KAKVSLITSD GIKFPQDGT L EFS DVTVDQT TGSITLRAIF PNP DHTLLPG MFVRRARLEEG LNP NAILVPQ QGVTRTPRGD ATVLVVGADD KVETRPIVAS QAIGDKWLVT EGLKAGDRVV ISGLQKVRPG VQVKAQEVTA DNNQQAASGA QPEQSKS
Source	Yeast
Target Names	acrA
Protein Names	Recommended name: Acriflavine resistance protein A
Expression Region	25-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.